



FELASA PANEL - SAMPLES					
IMMUNOCOMPETENT RAT					
Sample	ANNUAL	QUARTERLY	BASIC ANNUAL	BASIC QUARTERLY	
Hantaan virus	ELISA		ELISA		
Kilham rat virus	ELISA	ELISA	ELISA		ELISA
Lymphocytic choriomeningitis virus (LCMV)	ELISA				
Murine adenovirus F1 (MAD-1)	ELISA		ELISA		
Murine adenovirus K87 (MAD-2)	ELISA		ELISA		
Pneumonia virus of mice (PVM)	ELISA	ELISA	ELISA		ELISA
Rat Parvoviruses	ELISA	ELISA	ELISA		ELISA
Rat Minute Virus	ELISA	ELISA	ELISA		ELISA
Rat coronavirus/ Sialoadenitis virus (RCV/SDAV)	ELISA	ELISA	ELISA		ELISA
Reovirus type 3 (REO-3)	ELISA		ELISA		
Sendai virus	ELISA		ELISA		
Theiler's murine virus (TMEV)	ELISA	ELISA	ELISA		ELISA
Toolan's H-1 virus (H-1)	ELISA	ELISA	ELISA		ELISA
Cilia-Associated Respiratory Bacillus (CARB)	ELISA		ELISA		
Mycoplasma pulmonis (Mpul)	ELISA	ELISA	ELISA		ELISA
Clostridium piliforme – Tyzzer's (Cpil)	ELISA	ELISA	ELISA		ELISA
Streptobacillus moniliformis	ELISA		ELISA		
Helicobacter spp	PCR (SEQ)	PCR (SEQ)	PCR (SEQ)		PCR (SEQ)
Citrobacter rodentium	CULTURE	CULTURE	CULTURE		CULTURE
Corynebacterium kutscheri	CULTURE	CULTURE	CULTURE		CULTURE
Salmonella spp	CULTURE	CULTURE	CULTURE		CULTURE
Bordetella bronchiseptica	CULTURE	CULTURE	CULTURE		CULTURE
Pasteurella pneumotropica	CULTURE	CULTURE	CULTURE		CULTURE
Streptococcus β-haemolyticus (excl. group D)	CULTURE	CULTURE	CULTURE		CULTURE
Streptococcus pneumoniae	CULTURE	CULTURE	CULTURE		CULTURE
Opportunistic and commensal bacteria	CULTURE	CULTURE			
Pneumocystis carinii	ELISA		ELISA		
Encephalitozoon cuniculi (ECUN)	ELISA				
Ectoparasites	Tape	MICROSCOPY	MICROSCOPY		MICROSCOPY
Endoparasites - pinworms	Faeces	PCR	PCR		PCR

Sampling kit:

Dry Blood Spot Card	blood
Swab with medium	mouth and palate
Tube	faeces
Tape + slide	skin - hair

Storage conditions before shipment: 4°C