

	IMMUNOSUPRESSED MOUSE - LIVE ANIMAL PROFILES					
	FELASA PLUS		FELASA		BASIC FELASA	
	ANNUAL	QUARTERLY	ANNUAL	QUARTERLY	ANNUAL	QUARTERLY
Necropsy	●	●	●	●	●	●
HISTOPATHOLOGY (different organs)	●	●				
Ectromelia virus (ECTRO)	PCR		PCR		PCR	
Epizootic diarrhea of infant mice (EDIM)	RT-PCR	RT-PCR	RT-PCR	RT-PCR	RT-PCR	RT-PCR
Hantaan virus	RT-PCR		RT-PCR			
K virus (Kvir)	PCR		PCR			
Lymphocytic choriomeningitis virus (LCMV)	RT-PCR		RT-PCR		RT-PCR	
Murine adenovirus FI (MAD-1)	PCR		PCR		PCR	
Murine adenovirus K87 (MAD-2)	PCR		PCR		PCR	
Mouse cytomegalovirus (MCMV)	PCR		PCR			
Mouse hepatitis virus (MHV)	RT-PCR	RT-PCR	RT-PCR	RT-PCR	RT-PCR	RT-PCR
Murine parvovirus (NS1)	PCR	PCR	PCR	PCR	PCR	PCR
Mouse parvovirus (MPV 1 y 2)						
Minute virus of mice (MVM)	PCR	PCR	PCR	PCR	PCR	PCR
Murine norovirus (MNV)	RT-PCR	RT-PCR	RT-PCR	RT-PCR	RT-PCR	RT-PCR
Polyoma virus (Poly)	PCR		PCR			
Pneumonia virus of mice (PVM)	RT-PCR	RT-PCR	RT-PCR	RT-PCR	RT-PCR	
Reovirus type 3 (REO-3)	RT-PCR		RT-PCR		RT-PCR	
Sendai virus	RT-PCR	RT-PCR	RT-PCR	RT-PCR	RT-PCR	
Theiller's murine encephalomyelitis virus (TMEV)	RT-PCR	RT-PCR	RT-PCR	RT-PCR	RT-PCR	RT-PCR
Murine thymic virus (MTLV)	PCR		PCR			
Lactate dehydrogenase elevating virus (LDHV)	RT-PCR		RT-PCR			
Cillia-Associated Respiratory Bacillus (CARB)	PCR		PCR			
Mycoplasma pulmonis (Mpul)	PCR		PCR		PCR	
Clostridium piliforme – Tyzzer's (Cpil)	PCR/HISTOPATH	HISTOPATH	PCR		PCR	
Helicobacter spp	PCR (SEQ)	PCR (SEQ)	PCR (SEQ)	PCR (SEQ)	PCR (SEQ)	PCR (SEQ)
Bordetella bronchiseptica	CULTURE	CULTURE	CULTURE	CULTURE		
Citrobacter rodentium	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE
Corynebacterium kutscheri	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE
Pasteurella pneumotropica	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE
Salmonella spp	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE
Streptococcus β-haemolyticus (excl. group D)	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE
Streptococcus pneumoniae	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE	CULTURE
Streptobacillus moniliformis	PCR/HISTOPATH	HISTOPATH	PCR		PCR	
Opportunistic and commensal bacteria	CULTURE	CULTURE	CULTURE	CULTURE		
Pneumocystis murina	PCR/HISTOPATH	HISTOPATH	PCR			
Encephalitozoon cuniculi (ECUN)	PCR		PCR			
Ectoparasites	PARASITOLOGY	PARASITOLOGY	PARASITOLOGY	PARASITOLOGY	PARASITOLOGY	PARASITOLOGY
Endoparasites	PARASITOLOGY	PARASITOLOGY	PARASITOLOGY	PARASITOLOGY	PARASITOLOGY	PARASITOLOGY